

Glass DC Nebulizers

SeaSpray

2mL/min uptake	A13-07-USS2
400µL/min uptake	A13-07-USS04

Conikal

2mL/min uptake	A13-07-UC2
1mL/min uptake	A13-07-UC1

Slurry

4mL/min uptake	A13-07-US6
----------------	------------



DC Nebulizers for HF

DuraMist

1mL/min uptake	A13-07-DM1
----------------	------------



OpalMist

2mL/min uptake	A13-07-PFA2
----------------	-------------



VeeSpray

0.6 - 3mL/min uptake	A13-07-CV6
----------------------	------------



IsoMist XR

Programmable Temperature Cyclonic Spray Chamber

IsoMist XR Kit	KT-1137-XR
----------------	------------

Comprises:


Glass Twister Spray Chamber	21-809-3186
Assembly (ROL)	20-809-5176
PuraFlex Transfer Tube	21-809-4038
Metal Ball Joint Clip	JC-13
UniFit Right Angle Connector 2.0mm (PKT 3)	GAZ-04-2-90



BC Spray Chambers

Tracey BC

30mL Cyclonic	20-809-5316
* Direct replacement for 20-808-8882HE	



Tracey BC PEEK

30mL Cyclonic	20-809-4795
* Direct replacement for 20-809-2507	
* No internal resurfacing required	



Spray Chambers

Twister

50mL Cyclonic	20-809-9199HE
---------------	---------------



Tracey

50mL Cyclonic	20-808-8882HE
---------------	---------------



HydraMist Simultaneous Cold Vapor/ Pneumatic Nebulization Spray Chamber

Hydride Kit	KT-1157
-------------	---------



Spray Chambers for HF

Twister TFE Spray Chamber

50mL Cyclonic	20-809-3452
---------------	-------------



Tracey TFE Spray Chamber

50mL Cyclonic	20-809-2507
---------------	-------------



Torches

E-Torch

SVDV/VDV 30-808-4466

Comprises:

Body 31-808-4487
 Quartz Plasma Tube 31-808-4471
 Quartz Injector 1.8mm 31-808-3576
 Retainer 31-808-4589
 Ball Joint Adaptor 31-808-4565



NEW*

E-Torch

RV 30-808-4560

Comprises:

Body 31-808-4487
 Quartz Plasma Tube 31-808-4562
 Quartz Injector 1.8mm 31-808-3576
 Retainer 31-808-4589
 Ball Joint Adaptor 31-808-4565



NEW*


E-Torch Options:

Ceramic Plasma Tube for SVDV/VDV 31-808-4588
 Ceramic Plasma Tube for RV 31-808-4592
 Quartz Plasma Tube for SVDV/VDV 31-808-4471
 Quartz Plasma Tube for RV 31-808-4562
 O-ring kit 70-803-2921



E-Torch Tapered Injector Options

Quartz Injector 0.8mm 31-808-3738
 Quartz Injector 1.4mm 31-808-3583
 Quartz Injector 1.8mm 31-808-3576
 Quartz Injector 2.4mm 31-808-3584
 Alumina Injector 1.5mm 31-808-3725
 Alumina Injector 1.8mm 31-808-3585
 Alumina Injector 2.4mm 31-808-3722



Injectors for Agilent FDT

Quartz Injector 0.8mm 31-808-4357
 Quartz Injector 1.4mm 31-808-4354
 Quartz Injector 1.8mm 31-808-4355
 Quartz Injector 2.4mm 31-808-4356
 Alumina Injector 0.8mm 31-808-4530
 Alumina Injector 1.4mm 31-808-4532
 Alumina Injector 1.8mm 31-808-4358
 Alumina Injector 2.4mm 31-808-4534




Torches

(continued)

Quartz Torch

SVDV/VDV 31-808-3557
 RV 31-808-3556
 Retaining Ring 31-808-3563
 Pre-Optic Window Axial 70-803-1968
 Pre-Optic Window Radial 70-803-1969
 O-ring for Window Axial (PKT 5) 70-803-2794
 Pre-optics cone 70-803-2751



RF Coils

Silver 70-900-5100S
 Gold 70-900-5100G
 Copper 70-900-5100C
 RF Coil Installation Kit 70-900-5100T



RotaValve

RotaValve 70-803-1738



Valve Tubing Kit 70-803-1869



Sample Loop 1.0mL 70-803-1860



Nebulizer Connector UniFit - 1/4-28 Flangless 70-803-1888



Common Accessories

Elegra Argon Humidifier

Elegra Argon Humidifier 70-803-1268



Guardian In-Line Filters

Guardian DC in-line non-return gas filter 70-803-1942

Guardian in-line sample filter 70-803-1108


Eluo Adaptor for Inline Filter 70-803-1160



Trident Internal Standard Kits

Trident CT Standard Kit 60-703-1179

Standard Kit for HF 60-808-1150



Magnifier Inspection Tool 70-803-1923



Autosampler Probes

(continued)


For SPS 3/SPS 4

Guardian Probe 0.75mm

Probe Connecting Line (Red) 70-803-1957

PTFE Encapsulated Carbon Fibre Probe
0.75mm ID 70-803-0908

PTFE Encapsulated Carbon Fibre Probe
1.0mm ID 70-803-0853



Autosampler Probes (See page 22-23 for complete list of probes)

For Cetac ASX-200/500/800 Series

Guardian Probe 0.75mm

Probe Connecting Line (Red) 70-803-1803

Polyimide sheathed Autosampler Probe
0.5mm ID 60-808-1186L

PTFE Encapsulated Carbon Fibre Probe
0.5mm ID 70-803-0784

PTFE Encapsulated Carbon Fibre Probe
0.75mm ID 70-803-1443

PTFE Encapsulated Carbon Fibre Probe
1.0mm ID 70-803-0793



For Cetac ASX-7400, 7600 Series

Guardian Probe 0.75mm

Probe Connecting Line (Red) 70-803-2940

