

Glass DC Nebulizers

SeaSpray

400µL/min uptake

A23-1-USS04

MicroMist

400µL/min uptake

A23-1-UM04

Conikal

1mL/min uptake

A23-1-UC1



Quartz DC Nebulizers

SeaSpray

1mL/min uptake

A23-1-QSS1



DC Nebulizers for HF

DuraMist

400µL/min uptake

A23-1-DM04



OpalMist

400µL/min uptake

A23-1-PFA04



Spray Chambers

Twister

50mL Cyclonic

20-809-4417HE

50mL Cyclonic

20-809-4560

With Matrix Port



Tracey

50mL Cyclonic

20-809-4419HE



Spray Chambers

Cinnabar

20mL Cyclonic


20-809-4480HE



High Efficiency Sample Introduction System

HE-SIS for PerkinElmer NexION

KT-1184



Torches

D-Torch

30-808-3927



Comprises:

Base and Inner Tube

31-808-3928

Quartz Outer Tube

31-808-3942

Retaining Ring

31-808-2818

Injector Adaptor

31-808-3980

Injector Ferrule

31-808-2962

Fastener Assembly

31-808-4130

Viton O-rings (PKT 10)

70-V-011

PTFE Coated Viton O-rings
(PKT 10)

70-VT-012

Quartz Injectors

Tapered 1.0mm

31-808-3423

Tapered 1.8mm

31-808-4005

Tapered 2.0mm

31-808-4006

Tapered 2.5mm

31-808-3302



HF Resistant Injectors

Tapered Alumina 2.4mm

31-808-3205

Sapphire 2.4mm

31-808-3378



Cones

Sampler Cones

Nickel Sampler PE3011-Ni



Aluminium Sampler PE3011-Al

Platinum Sampler PE3013-Pt

Platinum Sampler (18mm) PE3013A-Pt

Platinum Sampler (18mm) PE3013C-Pt

Platinum Sampler Boron Free PE3013-Pt-BF



Skimmer Cones


Nickel Skimmer PE3012-Ni



Aluminium Skimmer PE3012-Al

Platinum Skimmer PE3014-Pt

Platinum Skimmer Boron Free PE3014-Pt-BF



Aluminium Hyper Skimmer PE3015-Al

Nickel Hyper Skimmer PE3015-Ni

Platinum Hyper Skimmer PE3015-Pt/Ni

O-ring for NexION Hyper Skimmer Cone (pack of 5) PE5115



Screw for NexION Hyper Skimmer Cone (pack of 2) PE5116



ConeGuard Thread Protectors

Sampler Cone Guard 70-803-1024
For NexION



Skimmer Cone Guard 70-803-1026
For NexION



Gasket for Sampler Cone (pack of 5) PE5111



Common Accessories

Elegra Argon Humidifier

Elegra Dual Argon Humidifier 70-803-1473



Guardian In-Line Filters

Guardian DC in-line non-return gas filter 70-803-1942

Guardian in-line sample filter 70-803-1108


Eluo Adaptor for Inline Filter 70-803-1160



Trident Internal Standard Kits

Trident CT Standard Kit 60-703-1179

Standard Kit for HF 60-808-1150



Magnifier Inspection Tool 70-803-1923



Autosampler Probes (See page 22-23 for complete list of probes)

For PerkinElmer S20 Series

Guardian Probe 0.75mm

Probe Connecting Line (Red) 70-803-2097



For Cetac ASX-200/500/800 Series

Guardian Probe 0.75mm

Probe Connecting Line (Red) 70-803-1803

Polyimide sheathed Autosampler

Probe 0.5mm ID 60-808-1186L

PTFE Encapsulated Carbon Fibre

Probe 0.5mm ID 70-803-0784

PTFE Encapsulated Carbon Fibre

Probe 0.75mm ID 70-803-1443

PTFE Encapsulated Carbon Fibre

Probe 1.0mm ID 70-803-0793

