

### Glass DC Nebulizers

**SeaSpray**

400µL/min uptake      A23-1-USS04

**MicroMist**

400µL/min uptake      A23-1-UM04

**Conikal**

1mL/min uptake      A23-1-UC1



### Quartz DC Nebulizers

**SeaSpray**

1mL/min uptake      A23-1-QSS1



### DC Nebulizers for HF

**DuraMist**

400µL/min uptake      A23-1-DM04



**OpalMist**

400µL/min uptake      A23-1-PFA04



### Spray Chambers

**Twister**

50mL Cyclonic      20-809-4417HE

50mL Cyclonic      20-809-4560

With Matrix Port



**Tracey**

50mL Cyclonic      20-809-4419HE



### Spray Chambers


**Cinnabar**

20mL Cyclonic      20-809-4480HE



**High Efficiency Sample Introduction System**

HE-SIS for PerkinElmer NexION      KT-1184



### Torches

**D-Torch**

30-808-3927

**Comprises:**

- Base and Inner Tube      31-808-3928
- Quartz Outer Tube      31-808-3942
- Retaining Ring      31-808-2818
- Injector Adaptor      31-808-3980
- Injector Ferrule      31-808-2962
- Fastener Assembly      31-808-4130
- Viton O-rings (PKT 10)      70-V-011
- PTFE Coated Viton O-rings (PKT 10)      70-VT-012



**Quartz Injectors**

- Tapered 1.0mm      31-808-3423
- Tapered 1.8mm      31-808-4005
- Tapered 2.0mm      31-808-4006
- Tapered 2.5mm      31-808-3302



**HF Resistant Injectors**

- Tapered Alumina 2.4mm      31-808-3205
- Sapphire 2.4mm      31-808-3378



**Cones**

**Sampler Cones**

Nickel Sampler PE3011-Ni



Platinum Sampler PE3013-Pt

Platinum Sampler (18mm) PE3013A-Pt

Platinum Sampler (12mm) PE3013C-Pt

Platinum Sampler PE3013-Pt

Boron Free PE3013-Pt-BF



**Skimmer Cones**

Nickel Skimmer PE4012-Ni

Platinum Skimmer PE4014-Pt

Nickel Hyper Skimmer PE4015-Ni

All Nickel Hyper Skimmer PE4016-Ni



**O-ring for NexION Hyper Skimmer Cone (pack of 5)**

PE5125



**ConeGuard Thread Protectors**

Sampler Cone Guard For NexION 70-803-1024



**Gasket for Sampler Cone (pack of 5)**

PE5111



**Common Accessories**

**Elegra Argon Humidifier**

Elegra Dual Argon Humidifier 70-803-1473



**Guardian In-Line Filters**

Guardian DC in-line non-return gas filter 70-803-1942

Guardian in-line sample filter 70-803-1108


Eluo Adaptor for Inline Filter 70-803-1160



**Trident Internal Standard Kits**

Trident CT Standard Kit 60-703-1179

Standard Kit for HF 60-808-1150



**Magnifier Inspection Tool 70-803-1923**



**Autosampler Probes** (See page 22-23 for complete list of probes)

**For PerkinElmer S20 Series**

Guardian Probe 0.75mm Probe Connecting Line (Red) 70-803-2097



**For Cetac ASX-200/500/800 Series**

Guardian Probe 0.75mm Probe Connecting Line (Red) 70-803-1803

Polyimide sheathed Autosampler Probe 0.5mm ID 60-808-1186L

PTFE Encapsulated Carbon Fibre Probe 0.5mm ID 70-803-0784

PTFE Encapsulated Carbon Fibre Probe 0.75mm ID 70-803-1443

PTFE Encapsulated Carbon Fibre Probe 1.0mm ID 70-803-0793

